

Savoir-Vivre

Curious how the child care bed
Savoir-Vivre 100 x 200 cm
(100 cm barred side height)
looks in your nursery?



Simply scan the QR code and
place the bed visually in your
nursery!



Strong as an animal through everyday life – our ergonomic beds offer your child the ideal conditions for a healthy sleep.

This model combines safety with comfort. The electric lying-height adjustment is especially useful when it comes to the daily care of handicapped people.



Electric height adjustment



Head- & footrest manual or electric adjustment (option)



EN 60601-2-52
DIN EN 50637

According to Medical Device Regulation
CE-Markering
Declaration of Conformity

Savoir-Vivre

Why a bed from SAVI?

- Great safety for your child
- Biologically neutral, pollutant-free (DIN EN 71-3)
- Stability and reliability
- Comfortable, back-friendly care height electric
- Many bed sizes
- Individual equipment options
- Optimal eye contact with your child
- Large, safe indoor area
- Colourful padding for maximum protection
- Short delivery times



Barred side height 125 cm,
e.g. in orange colour accent



Barred side height 60 cm,
e.g. in green colour accent



Barred side height 150 cm,
e.g. in purple colour accent



Barred side height 100 cm,
e.g. in blue colour accent



Barred side height 170 cm,
e.g. in blue-purple colour accent

Savoir-Vivre's individual strengths

- Barred side height 60 - 220 cm (option)
- Trendelenburg-/Antitrendelenburg-function (option)
- Drawer / storage (option)

All SAVI-locks
available for even
greater security

Technical data

Lying area	External dimension	Article-no.
100 x 200 cm	113 x 206 cm	
90 x 200 cm (size L)	103 x 206 cm	
90 x 170 cm (size M)	103 x 176 cm	
65 x 155 cm (size S)	78 x 161 cm	
Doors		
Four-winged folding sliding door		•
Four-winged folding sliding door on both sides	99.750	+
Folding door in the 270° version	99.704	+
Folding door in the 270° version on both sides	99.704BEID	+
Locking mechanism		
Child-proof SAVI-lock I		•
Child-proof SAVI-lock I metal	40100003	+
SAVI-lock II	99.033	+
SAVI-lock III	99.034	+
SAVI-lock IV	99.034A	+
Lying height (top of slatted bed base)		
40 - 80 cm electric height adjustment (without castors)		•
Barred side height		
	Overall height	
60/70 cm	110 cm	110.50.504.99101
100/110 cm	150 cm	150.50.504.99101
	for size 90 x 170 cm	
125/135 cm	175 cm	175.50.504.99101
	for size 90 x 170 cm	
150/160 cm	200 cm	200.50.504.99101
170/180 cm	220 cm	220.50.504.99101
Special sizes on request		
Slatted bed base		
Manual adjustment	99.030	•
Electric adjustment	99.200	+
Trendelenburg-/ Antitrendelenburg-funktion (tilting of the entire bed)	99.213	+
Running gear		
Four castor running gear with total locking mechanism (height: 17 cm)	99.887	+
Designer wooden castors (height: 17 cm)	99.886	+
Cut out on the front side (height: 8 cm)	99.830	n.a.c.
Kind of wood		
Pinewood		n.a.c.
Beechwood „natural oiled“ (further colours see p. 79)		•
Mattresses		
Polyfoam mattress	99.031.4	•
Cold foam mattress	99.031.5	+
Reha-star mattress plus (up to stage 2)	99.031.3.PLUS	+
Reha-star mattress II (up to stage 4)	99.031.3.ZWEI	+
Incontinence cover	99.031.3ZA50	+

• = standard o.r. = on request + = accessories n.a.c. = no additional charge

Technical data

Coloured features		
Natural	90.100	n.a.c.
Yellow	90.128	n.a.c.
Orange	90.204	n.a.c.
Blue-purple	90.405	n.a.c.
Purple	90.406	n.a.c.
Blue	90.515	n.a.c.
Green (for more colours see page 78)	90.618	•
Equipment		
Children's name on side of bed	Kindername	n.a.c.
Coloured sticks		n.a.c.
Perspex instead of sticks		o.r.
Tube feed-through	99.840	o.r.
External lock on side of bed	99.706	+
Raising aid	99.802	+
Sit up aid	99.820	+
Net-roof	99.800	+
Adaption of the bed for sloping ceiling		o.r.
Mattress holder	99.891	+
Protection for the bed base in case of kicks	99.892	+
Wood care set	99.880	+
Padding		
Motifpadding MALU, 50 cm	99.050MALU	+
Skai-padding 50 cm height	99.051	+
Skai-padding 100 cm height	99.051.100	+
Fixed Skai-padding	99.051.FEST	o.r.
Fixed Skai-padding with vinyl windows	99.051.FEST	o.r.
Edge padding	99.053	+

• = standard o.r. = on request + = accessories n.a.c. = no additional charge

Our Skai-padding is in many colours available!

You can find a pre-selection on our website at:

www.savi.de/en/bed-padding/



Savoir-Vivre

Key data of Savoir-Vivre

- Barred side height 60 - 170 cm, other heights on request
- Electric lying height 40 - 80 cm
- Four-winged folding sliding door
- Child-proof SAVI-lock
- Four-casters running gear (option)
- Perspex (option)

Heights

Barred side	Door	Total
• 60 cm	• 70 cm	• 110 cm
• 100 cm	• 110 cm	• 150 cm
• 125 cm	• 135 cm	• 175 cm
• 150 cm	• 160 cm	• 200 cm
• 170 cm	• 180 cm	• 220 cm



Total heights: 110 cm | 150 cm | 175 cm | 200 cm | 220 cm



Locking of doors



External lock on side of bed (option)

heavy, massive, good

Cut out



Height: 8 cm
(other heights are possible, e.g. for use of patient lifts)

great for the back



Examples



ID: 32835
Folding sliding door on both sides



ID: 34822
Coloured sticks



ID: 33458
White lacquered



ID: 34982
Completely perspex and drawers



ID: 34467
Foot cut-out (15 cm)
in the front frame

Savoir-Vivre

Optimum protection thanks to fixed Skai-padding

Your child is very active? With the firmly integrated Skai-padding (75 mm thickness), you offer your child optimum protection.

For more information on the different padding options, see page 79/80!



ID: 36814



ID: 32390
Stained white



ID: 36658



ID: 37345
Net-roof



ID: 33646

Our Skai-padding is in many colours available!
 You can find a pre-selection on our website at:
www.savi.de/en/bed-padding/



Optimal view thanks to integrated vinyl windows

The integrated vinyl windows fit flush with the fixed Skai-padding and thus offer a similar padding effect to the permanently integrated foam padding. The windows offer an optimal view. So you keep an eye on your child easily and your child feels visibly more comfortable.



ID: 34909



ID: 31485
Pine wood



ID: 23297



ID: 32813
Net-roof

Create your own



Create your children's bed for rehabilitation!

This is a list which you can use to create your children's bed for rehabilitations. We will take your wishes into consideration and make an offer about one of our bed's.

Just make a cross and decide how your bed should look like:

1. Lying area of the bed?

- 65 x 155 cm 70 x 140 cm 90 x 170 cm
 90 x 200 cm 100 x 200 cm other size: _____

2. Which barred side do you need? Height (cm): _____

3. Do you use a patient-lifter? Yes No

4. Which lying height (top of slatted base) do you need?

- 40 cm 50 cm 60 cm 70 cm other height: _____

5. Do you need an electric height adjustment? Yes No

6. Which slatted bed base do you need?

- no adjustment manual adjustment step-protection for the bed base
 electric adjustment Trendelenburg-function

7. What kind of doors do you want?

- 4-winged folding sliding door 270° door side rail "hide away" folding door
 lowerable side rail electric lowerable side rail

8. How many doors do you need?

- front side front and back sides head and footelement

9. Locking mechanism?

- SAVI-lock I SAVI-lock II SAVI-lock III SAVI-lock IV

10. Do you need perspex?

- Yes. Where?: _____ No

11. Kind of wood / wood-colour? Wood:

- beech pine
Wood-colour: natural basalt wenge barrique teak walnut
 white stained white lacquered

12. Colour feature?

- natural yellow orange red grey blue-purple purple
 blue green white pink other colour: _____
 coloured sticks children's name: _____

13. Running gear? castors without castors

14. Do you want drawer? Yes No

15. Do you want a mattress? Yes No

If yes, which one?

- polyfoam mattress interior spring mattress cold foam mattress
 Reha-Star II Reha-Star Plus incontinence cover

16. Do you need a padding? Yes No

If yes, which one?

- loose Skai-padding fixed Skai-padding edge-padding motifpadding MALU
height of the padding: _____ colour of the padding: _____

17. Do you need Vinyl-windows?

- Yes. Where? _____ No

18. Equipment?

- removable posts (Niklas) net-roof raising aid lockable handset box
 tube feed-through (gap) infusion holder drawers mattress holder
 touch-protection for height adjustment foot cut-out (e.g. for patient lift)
 central locking device for 4-roller chassis (Knut a+b only) access aid made of artificial leather
 supply opening (hole) with closing mechanism

Please use this statements to create an offer, where you choose a model, which fulfils the handicaps to the greatest possible extent.

Please use this statements to create an offer about a special model to fulfil the handicaps exactly.

place, date

stamp, sign

Very easy:

Just send us your order via email

or use our online-form at: www.savi.de/en/home/

SAVI GmbH Reha Solutions

Lise-Meitner-Str. 10 Fon: +49 (0) 25 64 - 93 09-70

D-48691 Vreden Fax: +49 (0) 25 64 - 93 09-60