

Lola



## Small heroes very big – Lola is the ideal bed for toddlers.

This bed is ideal for small children because it has a fixed, high basic lying height, making nursing more comfortable. It also relieves pressure on the back during care and ensures your child has a restful night's sleep.



More Info about  
Lola:  
[www.savi.de/en/care-beds-for-children/](http://www.savi.de/en/care-beds-for-children/)

EN 60601-2-52  
DIN EN 50637

According to Medical Device Regulation  
CE-Markering  
Declaration of Conformity

## Why a bed from SAVI?

- Great safety for your child
- Biologically neutral, pollutant-free (DIN EN 71-3)
- Stability and reliability
- Comfortable, back-friendly care height electric
- Many bed sizes
- Individual equipment options
- Optimal eye contact with your child
- Large, safe indoor area
- Colourful padding for maximum protection
- Short delivery times

## Lola's individual strengths

- Good care height/slatted bed base 85 cm
- Drawer (option)



Choice of colours (see on page 78)

door open  
→



Barred side height 65 cm  
e.g. In green colour accent



Barred side height 65 cm and drawers  
e.g. In green colour accent

## Technical data

Lying area	External dimension	Article-no.	
100 x 200 cm	113 x 206 cm		
90 x 200 cm	103 x 206 cm		
<b>Doors</b>			
Four-winged folding sliding door			•
Four-winged folding sliding door on both sides	99.750		+
Folding door in the 270° version	99.704		+
Folding door in the 270° version on both sides	99.704BEID		+
<b>Locking mechanism</b>			
Child-proof SAVI-lock I			•
Child-proof SAVI-lock I metal	40100003		+
SAVI-lock II	99.033		+
SAVI-lock III	99.034		+
SAVI-lock IV	99.034A		+
<b>Lying height (top of slatted bed base)</b>			
85 cm (without castors)			•
<b>Barred side height</b>	<b>Overall height</b>		
65 cm	150 cm	150.85.002	•
90 cm	175 cm	175.85.002	+
Special sizes on request			+
<b>Slatted bed base</b>			
Manual adjustment	99.030		+
Electric adjustment	99.200		+
<b>Running gear</b>			
Four castor running gear with total locking mechanism	99.887		o.r.
Designer wooden castors	99.886		o.r.
<b>Kind of wood</b>			
Pinewood			n.a.c.
Beechwood „natural oiled“ (further colours see p. 79)			•
<b>Mattresses</b>			
Polyfoam mattress	99.031.4		•
Cold foam mattress	99.031.5		+
Reha-Star Plus (Stadium 2)	99.031.3 PLUS		+
Reha-Star II (Stadium4)	99.031.3.ZWEI		+
Incontinence cover	99.031.3.ZA50		+

• = standard    o.r. = on request    + = accessories    n.a.c. = no additional charge

## Technical data

<b>Coloured features</b>		
Natural	90.100	n.a.c.
Yellow	90.128	•
Orange	90.204	n.a.c.
Blue-purple	90.405	n.a.c.
Purple	90.406	n.a.c.
Blue	90.515	n.a.c.
Green (for more colours see page 78)	90.618	n.a.c.
<b>Equipment</b>		
Children's name on side of bed	Kindename	n.a.c.
Coloured sticks to go to the coloured feature		n.a.c.
Perspex instead of sticks		o.r.
External lock on side of bed	99.706	+
Tube feed-through	99.840	o.r.
Raising aid	99.802	+
Sit up aid	99.820	+
Adaption of the bed for sloping ceiling		o.r.
Drawer (2 pieces)	99.041	+
Mattress holder	99.891	+
Protection for the bed base in case of kicks	99.892	+
Wood care set	99.880	+
<b>Padding</b>		
Motifpadding MALU, 50 cm	99.050MALU	+
Skai-padding 50 cm height	99.051	+
Fixed Skai-padding	99.051.FEST	o.r.
Fixed Skai-padding with vinyl windows	99.051.FEST	o.r.
Edge padding	99.053	+

• = standard    o.r. = on request    + = accessories    n.a.c. = no additional charge



# Lola

## Key data of Lola

- Barred side height of 65 - 90 cm, other heights on request
- Fixed lying height 85 cm
- Child-proof SAVI-lock
- Drawer (option)
- Perspex (option)



ID: 17042  
Fixed Skai-padding in white



ID: 33029  
270° door with Lexan, the rest closed with wood instead of sticks



ID: 28513  
Lying area 130 x 160 cm, barred side height 45 cm



ID: 26380  
Four lowerable side rails



ID: 31388  
Wood colour 'Wenge',  
bed gallow

"I like to dream  
about my friends  
and the beach."



ID: 19540  
Supply opening (hole)  
with closing mechanism



ID: 36779  
Drawer, stained white



ID: 31727  
Perspex allover

# Create your own



## Create your children's bed for rehabilitation!

This is a list which you can use to create your children's bed for rehabilitations. We will take your wishes into consideration and make an offer about one of our bed's.

Just make a cross and decide how your bed should look like:

### 1. Lying area of the bed?

- 65 x 155 cm    70 x 140 cm    90 x 170 cm  
 90 x 200 cm    100 x 200 cm   other size: \_\_\_\_\_

### 2. Which barred side do you need? Height (cm): \_\_\_\_\_

### 3. Do you use a patient-lifter? Yes No

### 4. Which lying height (top of slatted base) do you need?

- 40 cm    50 cm    60 cm    70 cm   other height: \_\_\_\_\_

### 5. Do you need an electric height adjustment? Yes No

### 6. Which slatted bed base do you need?

- no adjustment    manual adjustment    step-protection for the bed base  
 electric adjustment    Trendelenburg-function

### 7. What kind of doors do you want?

- 4-winged folding sliding door    270° door    side rail "hide away"    folding door  
 lowerable side rail    electric lowerable side rail

### 8. How many doors do you need?

- front side    front and back sides    head and footelement

### 9. Locking mechanism?

- SAVI-lock I    SAVI-lock II    SAVI-lock III    SAVI-lock IV

### 10. Do you need perspex?

- Yes. Where?: \_\_\_\_\_    No

### 11. Kind of wood / wood-colour? Wood:

- beech    pine  
Wood-colour:  natural    basalt    wenge    barrique    teak    walnut  
 white stained    white lacquered

### 12. Colour feature?

- natural    yellow    orange    red    grey    blue-purple    purple  
 blue    green    white    pink    other colour: \_\_\_\_\_  
 coloured sticks    children's name: \_\_\_\_\_

### 13. Running gear? castors   without castors

### 14. Do you want drawer? Yes No

### 15. Do you want a mattress? Yes No

If yes, which one?

- polyfoam mattress    interior spring mattress    cold foam mattress  
 Reha-Star II    Reha-Star Plus    incontinence cover

### 16. Do you need a padding? Yes No

If yes, which one?

- loose Skai-padding    fixed Skai-padding    edge-padding    motifpadding MALU  
height of the padding: \_\_\_\_\_ colour of the padding: \_\_\_\_\_

### 17. Do you need Vinyl-windows?

- Yes. Where? \_\_\_\_\_    No

### 18. Equipment?

- removable posts (Niklas)    net-roof    raising aid    lockable handset box  
 tube feed-through (gap)    infusion holder    drawers    mattress holder  
 touch-protection for height adjustment    foot cut-out (e.g. for patient lift)  
 central locking device for 4-roller chassis (Knut a+b only)    access aid made of artificial leather  
 supply opening (hole) with closing mechanism

Please use this statements to create an offer, where you choose a model, which fulfils the handicaps to the greatest possible extent.

Please use this statements to create an offer about a special model to fulfil the handicaps exactly.

place, date

stamp, sign

**Very easy:**

Just send us your order via email

or use our online-form at: [www.savi.de/en/home/](http://www.savi.de/en/home/)

SAVI GmbH Reha Solutions

Lise-Meitner-Str. 10   Fon: +49 (0) 25 64 - 93 09-70

D-48691 Vreden   Fax: +49 (0) 25 64 - 93 09-60