

Jonas



On a dream journey with Jonas – our handmade beds offer your child an unmistakable sleeping comfort and particularly beautiful dreams.

This bed offers a high degree of safety. Guaranteed by the high barred sides and our proven SAVI-lock.



EN 60601-2-52
DIN EN 50637

According to Medical Device Regulation
CE-Markering
Declaration of Conformity

Why a bed from SAVI?

- Great safety for your child
- Biologically neutral, pollutant-free (DIN EN 71-3)
- Stability and reliability
- Comfortable, back-friendly care height electric
- Many bed sizes
- Individual equipment options
- Optimal eye contact with your child
- Large, safe indoor area
- Colourful padding for maximum protection
- Short delivery times



ID: 11071
Barred side height 110 cm,
e.g. with drawers,
in blue colour accent

Jonas' individual strengths

- Cost-effective, as no height adjustment is required
- Barred side 60 - 220 cm (option)
- Drawer/storage (option)



Choice of colours (see on page 78)



Barred side height 100 cm,
e.g. in blue colour accent

Technical data

Lying area	External dimension	Article no.
100 x 200 cm	113 x 206 cm	
90 x 200 cm	103 x 206 cm	
Doors		
Four-winged folding sliding door		•
Four-winged folding sliding door on both sides	99.750	+
Folding door in the 270° version	99.704	+
Folding door in the 270° version on both sides	99.704BEID	+
Locking mechanism		
Child-proof SAVI-lock I		•
Child-proof SAVI-lock I metal	40100003	+
SAVI-lock II	99.033	+
SAVI-lock III	99.034	+
SAVI-lock IV	99.034A	+
Lying height (top of slatted bed base)		
50 cm (without castors)		•
Barred side height		Overall height
100 cm	150 cm	150.50.004
125 cm	175 cm	175.50.004
Special sizes on request		
Slatted bed base		
Manual adjustment	99.030	+
Electric adjustment (no height adjustment)	99.200	+
Rolling gear		
Four castor running gear with total locking mechanism (height: 17 cm)	99.887	o.r.
Designer wooden castors (height: 17 cm)	99.886	o.r.
Kind of wood		
Pinewood		n.a.c.
Beechwood „natural oiled“ (further colours see p. 79)		•
Mattresses		
Polyfoam mattress	99.031.4	•
Cold foam mattress	99.031.5	+
Reha-star mattress plus (up to stage 2)	99.031.3.PLUS	+
Reha-star mattress II (up to stage 4)	99.031.3.ZWEI	+
Incontinence cover	99.031.3ZA50	+

• = standard o.r. = on request + = accessories n.a.c. = no additional charge



Technical data

Coloured features		
Natural	90.100	n.a.c.
Yellow	90.128	n.a.c.
Orange	90.204	n.a.c.
Blue-purple	90.405	n.a.c.
Purple	90.406	n.a.c.
Blue	90.515	•
Green (for more colours see page 78)	90.618	n.a.c.
Equipment		
Children's name on side of bed	Kindename	n.a.c.
Coloured sticks to go to the coloured feature		n.a.c.
Perspex instead of sticks		o.r.
External lock on side of bed	99.706	+
Tube feed-through	99.840	o.r.
Raising aid	99.802	+
Adaption of the bed for sloping ceiling		o.r.
2 drawers	99.040	+
Mattress holder	99.891	+
Protection for the bed base in case of kicks	99.892	+
Wood care set	99.880	+
Padding		
Motifpadding MALU, 50 cm	99.050MALU	+
Skai-padding 50 cm height	99.051	+
Skai-padding 100 cm height	99.051.100	+
Fixed Skai-padding	99.051.FEST	o.r.
Fixed Skai-padding with vinyl windows	99.051.FEST	o.r.
Edge padding	99.053	+

• = standard o.r. = on request + = accessories n.a.c. = no additional charge

Our Skai-padding is in many colours available!

You can find a pre-selection on our website at:

www.savi.de/en/bed-padding/



Jonas

Key data of Jonas

- Barred side height 100 - 125 cm, or even higher (i.e. 220 cm) available
- Fixed lying height 50 cm (entrance height)
- Child-proof SAVI-lock
- Drawers (option)
- Perspex (option)



ID: 35578
270° doors, fixed padding,
sticks, castors and Lexan



Many wood colours possible (see p. 79)



ID: 32416
Drawers, castors



ID: 32149
Wood colour 'Walnut'



ID: 37134
Perspex

Examples



ID: 34012
Lexan equipment
(with air holes)
and net-roof



ID: 35831
Fixed Skai padding with
vinyl windows and child's name
on the side of the bed

"I don't put my head
in the sand -
I put it on a
comfortable pillow."



ID: 30590
Entry aid made from
imitation leather

Create your own



Create your children's bed for rehabilitation!

This is a list which you can use to create your children's bed for rehabilitations. We will take your wishes into consideration and make an offer about one of our bed's.

Just make a cross and decide how your bed should look like:

1. Lying area of the bed?

- 65 x 155 cm 70 x 140 cm 90 x 170 cm
 90 x 200 cm 100 x 200 cm other size: _____

2. Which barred side do you need? Height (cm): _____

3. Do you use a patient-lifter? Yes No

4. Which lying height (top of slatted base) do you need?

- 40 cm 50 cm 60 cm 70 cm other height: _____

5. Do you need an electric height adjustment? Yes No

6. Which slatted bed base do you need?

- no adjustment manual adjustment step-protection for the bed base
 electric adjustment Trendelenburg-function

7. What kind of doors do you want?

- 4-winged folding sliding door 270° door side rail "hide away" folding door
 lowerable side rail electric lowerable side rail

8. How many doors do you need?

- front side front and back sides head and footelement

9. Locking mechanism?

- SAVI-lock I SAVI-lock II SAVI-lock III SAVI-lock IV

10. Do you need perspex?

- Yes. Where?: _____ No

11. Kind of wood / wood-colour? Wood:

- beech pine
Wood-colour: natural basalt wenge barrique teak walnut
 white stained white lacquered

12. Colour feature?

- natural yellow orange red grey blue-purple purple
 blue green white pink other colour: _____
 coloured sticks children's name: _____

13. Running gear? castors without castors

14. Do you want drawer? Yes No

15. Do you want a mattress? Yes No

If yes, which one?

- polyfoam mattress interior spring mattress cold foam mattress
 Reha-Star II Reha-Star Plus incontinence cover

16. Do you need a padding? Yes No

If yes, which one?

- loose Skai-padding fixed Skai-padding edge-padding motifpadding MALU
height of the padding: _____ colour of the padding: _____

17. Do you need Vinyl-windows?

- Yes. Where? _____ No

18. Equipment?

- removable posts (Niklas) net-roof raising aid lockable handset box
 tube feed-through (gap) infusion holder drawers mattress holder
 touch-protection for height adjustment foot cut-out (e.g. for patient lift)
 central locking device for 4-roller chassis (Knut a+b only) access aid made of artificial leather
 supply opening (hole) with closing mechanism

Please use this statements to create an offer, where you choose a model, which fulfils the handicaps to the greatest possible extent.

Please use this statements to create an offer about a special model to fulfil the handicaps exactly.

place, date

stamp, sign

Very easy:

Just send us your order via email

or use our online-form at: www.savi.de/en/home/

SAVI GmbH Reha Solutions

Lise-Meitner-Str. 10 Fon: +49 (0) 25 64 - 93 09-70

D-48691 Vreden Fax: +49 (0) 25 64 - 93 09-60